

# DRUGWIPE II<sup>®</sup>

detects drugs-of-abuse.



**Detection of:**

- ▶ Cannabis
- ▶ Opiates
- ▶ Cocaine
- ▶ Amphetamines, Methamphetamines, Ecstasy
- ▶ Benzodiazepines

# Rapid test for the detection of drugs-of-abuse on surfaces, in sweat or in saliva

## Characteristics:

### ● Detection of:

- ▶ Cannabis
- ▶ Opiates
- ▶ Cocaine
- ▶ Amphetamines, Methamphetamines, Ecstasy
- ▶ Benzodiazepines

Further tests in preparation

### ● Lower detection limits:

- ▶ surfaces: **2–50 ng**
- ▶ saliva / sweat: **5–300 ng/ml**

### ● Ease of use

- ▶ collection of sample by wiping
- ▶ inserting wiping part into test device
- ▶ automatic transfer of drug sample onto test device
- ▶ dipping into water
- ▶ reading of the result after approx. 3–10 min

### ● Maximum reliability

- ▶ very small cross reactivity due to highly specific reagents (antibodies)
- ▶ internal control line

### ● Quick to use

- ▶ examination time: approx. 3–10 min

### ● Wide range of application

- ▶ traffic control
- ▶ recognition of drug abuse
- ▶ custom search operations
- ▶ illegal drug searches

The DRUGWIPE II® rapid tests are designed for on-the-spot use during drug control operations, supporting a suspicion of drug abuse, drug smuggling or dealing. A final confirmation requires appropriate laboratory tests.



**SECURETEC**  
DETEKTIONS-SYSTEME AG

Rosenheimer Landstr. 129 · D-85521 Ottobrunn  
Phone: ++49(0)89/62 78 96-51 · Telefax: 62 78 96-52  
E-Mail: [info@securetec.net](mailto:info@securetec.net)  
<http://www.securetec.net>

Securetec develops and sells test systems and provides professional service and advice for the implementation of test devices and your anti drug policy.